

X-RAY PROTECTIVE APRON

Maximum Protection, Maximum Flexibility


*Light weight for comfort
with good design*










X-RAY PROTECTIVE APRON

- 0.25mm Pb Equivalent
- 0.5 mm Pb Equivalent
- SIZE : Large : 44 x 24 inches (110cm x 60cm)
Medium : 40 x 24 inches (100cm x 60cm)
Small : 36 x 24 inches (90cm x 60cm)
- COLOR : Pink, Blue, Dark Blue
- Any Sizes are Available upon Request.





X-RAY Protective Aprons		
Buckle Type	Velcro Type	
		
Coat Type	Vest & Skirt Type	
		
<ul style="list-style-type: none"> - Pb Equivalent : 0.25 mmpb, 0.5mmpb - Colours : Dark Blue, Blue, Pink - Size : Small(S), Medium(M), Large(L) 		
Thyroid collar	Gonad Shields	Ovary Shields
		
<ul style="list-style-type: none"> - Pb Equivalent : 0.5mmpb - Colour : Dark Blue, Blue, Pink 		