



High Frequency Portable X-Ray 9020H.F



9020

High Frequency Portable Generator

Specifications

- Ultra light weight-13lbs
 - 40-90kV step /20mA, 18mA, 16mA
 - 1.60 kW Max Output
 - 100 kHz h.f switching inverter
 - 2 Year Warranty
-
- Dual laser pointers (26") & measuring tape
 - Built in collimator
 - Reversible digital display of kV & mAs
 - Soft touch user panel
 - Overload protection for kV,mA time and temp.
 - Electronic voltage compensation
 - Self manageable APR functions with 5 memory keys
 - Foam grips on each side of handle for easy handling
 - Laser pointers and collimator light activate with a touch of the exposure switch
 - Carrying case for easy transport
 - Ideal for equine applications
-

Authorized Dealer

JPI Healthcare Solutions, Inc.
52 Newtown Plaza • Plainview , NY 11803
Toll Free : 888-574-2001
Tel: 516-513-1330
Fax: 516-513-1331
www.jphealthcare.com