

# examvue CR PRO

for Podiatry



## Specifications:

- Cassette Sizes: 14x17", 10x12", 8x10"  
*\*System Package includes two: 14x17" Cassette/IPs*
- System Set-up Options: Tabletop • Vertical wall mount
- Resolution: 100um/200um
- Data Capture: 16-bits per pixel, 65,535 gray scale
- Dimensions (HxWxD): 5.3" x 17.5" x 31.3" (134 x 444 x 795 mm)
- Weight: 42.4lbs (19.5kg)



## Acquire Software:

- DICOM Send, DICOM Print, CD/DVD Burning
- Multi Frequency Image Processing
- Zoom, Cropping, Mark, Window Level
- User Defined LUT
- Saving of measurements & annotations
- Embedded viewer



(computer mount not included)

## Viewing License:

- Distance Comparison
- Image Stitching
- Cobb Angle Measurements
- PACS connectivity for Easy Image Comparison
- Save Measurements and Annotations
- Prints on DICOM and Regular Printers
- DICOM Query/Retrieve /Send SCU/SCP
- DICOM Grayscale Print SCU
- View via Network



## Regulatory approvals:

- FDA 510(k) Premarket Notification
- CE (0120)
- UL, cUL, FCC
- Health Canada Certification
- Not approved for mammography



Authorized Dealer

JPI Healthcare Solutions Inc.  
52 Newtown Plaza Plainview NY 11803  
Toll Free: 888-574-2001  
Tel: 516-513-1330  
Fax: 516-513-1331  
www.jpihealthcare.com