



Radiography System



DynaVue



Ultrasound



X-Ray Grid



Software



Detector



X-Ray Accessories




Diagnostics

2-In-1 Digital X-ray System

DynaVue is a robust, user-friendly digital radiography system designed to produce great X-ray images of companion animals. In development of DynaVue, the JPI team worked side-by-side with industry-leading veterinarians to optimize image quality and diagnostic results. DynaVue is powered by ExamVue software, which was developed specifically with veterinarians and their workflows in mind.



for Veterinary



Radiography System

Superior Digital Imaging & Increased Practice Efficiency

- 4-way floating table for X-ray positioning with minimal disturbance to your patient
- Easy-to-use operating console with large touchscreen monitor installed on tube stand
- High-sensitivity Cesium Iodide (CsI) detector technology
- Low dose, high frequency generator and remote capture feature
- Compact unit with wheels for convenient positioning for easy cleaning
- 220V, 15A power requirements for easy installation in any location




Added Value of Fluoroscopy

With the push of a button, this space and timesaving system can go from digital X-ray to fluoroscopy mode mid-exam, without any adjustment to your patient. Applications beyond digital X-ray are plentiful with DynaVue's fluoroscopy capabilities. Some of these revenue-generating procedures include:

- Minimally invasive fracture repair
- Intraluminal stenting for tracheal collapse treatment
- Treating urethral and colonic obstruction
- Vascular procedures
- Contrast studies of the gastrointestinal tract

Fluoroscopy Sample Image



Radiography Sample Images

Canine

Canine

ExamVue software includes:

- Patient registration program at front desk
- Imaging and annotation tools
- Ability to email images or send by DICOM
- Vertebral Heart Score
- Powerful user features such as: Image Import/Export, Image Comparison, Multi-view, Image Expansion, Full Screen, Reject/Accept, Image Capture
- DICOM cine-loop display to view fluoroscopy video inside your radiographic studies

Best Warranty & Lowest Cost of Ownership

Our DynaVue system includes:

- Detector
- Generator
- 4-way floating tabletop
- Computer with touchscreen monitor
- ExamVue software

JPI provides an unrivaled warranty and coverage program. Please contact us for more details.