



Detector

TruePortable



Ultrasound



X-Ray Grid



Software



Radiography
System



X-Ray
Accessories



Diagnostics

Lightweight, Easy-to-Use, Portable X-ray System.

Lighter, more maneuverable, and better protected than the competition, True Portable was designed specifically with portable provider technologists in mind. The True Portable can be moved in any direction for easy positioning and quick maneuvering in tight areas.



Detector

System advantages include:

- 8 wheels
 - 4 free-moving wheels which allow indoor and bedside movement in any direction
 - 2 large wheels for travel and climbing stairs
 - 2 wheels on the handrails for loading into a vehicle
- Powerful 120 kVp X-ray source to ensure great images every time
- Arrives ready-to-use from JPI. Unfold and be ready to take an X-ray immediately
- High sensitivity Csl technology provides superior quality images
- Send to multiple PACS automatically



High sensitivity Csl technology provides superior quality images



Unfold and be ready to take an X-ray immediately

True Portable combines portability, durability, power, and advanced image processing to provide exceptional images. We built True Portable with protection of your investment in mind. Both the flat panel detector and laptop have their own fully protected encasements, which come standard in this rugged portable package. Great protection will allow your team to go about their day knowing that your most valuable components are well-protected no matter where you take them.



2 wheels on the handrails for loading into a vehicle

JPI's unrivaled Coverage Program provides the lowest 5-years cost of ownership by combining:

- 5 years of Detector Drop Coverage
- 5-year Detector Warranty
- 5 years of Remote Software Support
- 3-year Warranty on Computer
- 2-year Warranty on Generator
- 1-year Warranty on Stand