



## Our Solutions

With superior image quality and intuitive imaging software, the ExamVue digital imaging system delivers fast image acquisition, easy workflow management, and lower dose for your patients.

The ExamVue system comes equipped with a strong suite of chiropractic tools, and is available in any configuration for your practice needs.

## JPI ExamVue DX Package



DR17T  
17"x17" Tethered



PC Workstation



LAN PACS



X-ray Grid



ExamVue Acquisition  
Software



DICOM  
Modality Worklist

## Optional Upgrade



DR14W  
14"x17" Wireless



Laptop  
Workstation



Cloud PACS

## 5-5-5 Program

JPI is not your average imaging solutions provider. Reach out to us! Our Account Managers are happy to discuss the imaging solutions and coverage options that are best for you.

- ▶ **5 years DROP COVERAGE PROGRAM**
- ▶ **5 years REMOTE SOFTWARE SUPPORT**
- ▶ **5 years WARRANTY**





## ExamVue DX Flat Panel Detectors

- ▶ Works with all x-ray generators, no matter the age
- ▶ Installs directly into your existing bucky stand
- ▶ Sharp, easily read images
- ▶ Cesium Direct-Deposit Technology
  - Lower patient dose
  - Great images with larger patients
  - No compromise on image quality
- ▶ Automatic Exposure Detection
  - No wiring into the generator needed



### DR17T


- ▶ Never leave your stand
- ▶ No worry, no fuss
- ▶ Reliable consistency

### DR14W

- ▶ Use like a cassette
- ▶ Long battery life
- ▶ Fast data transmission
- ▶ Wifi Connection

## ExamVue Acquisition Software

With state of the art and intuitive imaging tools, our ExamVue software produces clear, concise, and quality images with even the oldest generators and larger patients. ExamVue PACS comes standard with chiropractic diagnostic tools and can be expanded even further with the JPI Cloud PACS. Reliable and convenient, the ExamVue software seamlessly integrates into any professional practice.




## ExamVue PACS Software

### Advanced Chiropractic Tool-Set



Cervical Curve



Gonstead Toolset

#### Chiropractic Tool-Set

- Cervical Curve
- Lumbar Curve
- Cobb Angle
- Gonstead Toolset
- Straight line measurement
- Norberg Angle
- Ilium Analysis
- Image Stitching
- Angle Measurement
- Spine Labeling
- Baseline Measurement
- Draw Line/Arrow

### Additional Tools in Cloud PACS

Horizontal/Vertical Distance  
Distance from Vertical Line  
Distance from Horizontal Line  
Spinal Body Canal Ratio Line  
Horizontal Deflection Line  
Angles Between lines  
Atlas Plane Line  
George's Line  
Ray Line

Center Mass  
Antero & Retrolisthesis  
Logan Basic Marking  
Full Lateral Spine Marking  
Gonstead Pelvic Analysis  
Extended Cobb's Angle  
Vertebral Angle  
Multiple Angles  
Atlas Rotation

**User-Friendly** Design

**Simple & Fast** Workflow

**Convenient Diagnosis** Function

#### Contact us:

3555 Veterans Memorial Hwy, Unit D  
Ronkonkoma, NY 11779  
sales@jpihealthcare.com  
(516) 513-1330 Ext. 3



www.jpihealthcare.com