

EV Glassless Substrate 10" x 12" Flat Panel Detectors



ExamVue Glassless Flat Panel Detectors

Panel Detector Features



Works with All X-ray Generators



Automatic Exposure Detection



Sharp and Easily Read Images



Lower Patient Dose with CSI Tech



Great Images with Larger Patients




Installs Directly into Bucky Stand



No Wiring into the Generator Needed

ExamVue Duo Acquisition Software

With state of the art and intuitive imaging tools, our ExamVue Duo software produces clear, concise, and quality images with even the oldest generators. ExamVue PACS comes standard with diagnostic tools and can be expanded even further with the JPI Cloud PACS. Reliable and convenient, the ExamVue Duo software seamlessly integrates into any professional practice.



Contact us:

3555 Veterans Memorial Hwy, Unit D
Ronkonkoma, NY 11779
sales@jpihealthcare.com
(516) 513-1330 Ext. 3



www.jpihealthcare.com

TECHNICAL
SPECIFICATIONS



ExamVue Detectors
Sensor
Scintillator
Active Area
Pixel Array
Pixel Pitch
Image Quality
Spatial Resolution
A/D Conversion
Full Image Time
MTF
DQE
Communication
Communication Interface
Exposure Control
Mechanical
Dimension
Weight
Environmental
Operating Temperature
Operating Humidity
Storage Temp
Storage Humidity
IP Rating
Power
Max Power (Standby/Operating)
Power Supply
Frequency
Battery Autonomy

DR10FW
Csl
315.2 mm x 250.2 mm
3152 x 2502
100 μm
4.5 sec
2.4G/5GHz, 300 Mbps
362.1 x 269.2 x 15mm
2 kg
10°C - 35°C
5% - 90% RH
-20°C - 55°C
5% - 95% RH
IP67
18W