



Our Solutions

With superior image quality and intuitive imaging software, the ExamVue Duo digital imaging system delivers fast image acquisition, easy workflow management, and lower dose for your patients.

The ExamVue Duo system comes equipped with a strong suite of advanced tools, and is available in any configuration for your practice needs.

JPI ExamVue Duo Package



ExamVue
14"x17" Wireless



PC Workstation



LAN PACS



X-ray Grid



ExamVue Duo
Acquisition Software



DICOM
Modality Workflow

Optional Upgrade



ExamVue
17"x17" Tethered



Laptop
Workstation



Cloud PACS



CubeX 28
X-ray Generator



5-5-5 Warranty Program

JPI is not your average imaging solutions provider. Reach out to us! Our Account Managers are happy to discuss the imaging solutions and coverage options that are best for you.


- ▶ Up to 5 years of **HARDWARE WARRANTY**
- ▶ Up to 5 years of **DROP COVERAGE PROGRAM**
- ▶ Up to 5 years of **REMOTE SOFTWARE SUPPORT**





ExamVue Duo Flat Panel Detectors

- ▶ Works with all x-ray generators, no matter the age
- ▶ Installs directly into your existing bucky stand
- ▶ Sharp, easily read images
- ▶ Cesium Direct-Deposit Technology
 - Lower patient dose
 - Great images with larger patients
 - No compromise on image quality
- ▶ Automatic Exposure Detection
 - No wiring into the generator needed



17" x 17" Tethered


- ▶ Never leave your stand
- ▶ No worry, no fuss
- ▶ Reliable consistency

14" x 17" Wireless

- ▶ Use like a cassette
- ▶ Long battery life
- ▶ Fast data transmission
- ▶ Wifi Connection

ExamVue Duo Acquisition Software

With state of the art and intuitive imaging tools, our **ExamVue Duo** software produces clear, concise, and quality images with even the oldest generators and largest patients. Reliable and convenient, the **ExamVue Duo** software seamlessly integrates into any professional medical practice.



ExamVue Duo comes standard with medical diagnostic tools and can be expanded even further with the JPI Cloud PACS. The JPI Cloud PACS expands your access to your images to anywhere you can connect to the internet and allows for easy sharing with doctors and patients throughout your practice or anywhere in the world.

ExamVue Duo Advanced Tool-set



- Center Point
- Orthogonal Line
- Torg's Ratio
- Pelvic Obliquit
- Baseline Measurement
- Jackson Angle
- Ishihara Index
- George's Line
- Cervical Curve
- Lumbar Curve
- Ilium Analysis
- Spine Analysis
- Spine Label

Additional Tools in Cloud PACS

- | | |
|-------------------------------|----------------------------|
| Linear/Angle Measure | Atlas Plane Line |
| Circle/Polygon Area | George's Line |
| Text Annotation | Ray Line |
| Arrow/Line Mark | Center Mass |
| Horizontal/Vertical Distance | Antero & Retrolisthesis |
| Distance from Vertical Line | Full Lateral Spine Marking |
| Distance from Horizontal Line | Extended Cobb's Angle |
| Spinal Body Canal Ratio Line | Vertebral Angle |
| Horizontal Deflection Line | Multiple Angles |
| Angles Between lines | Atlas Rotation |

User-Friendly Design

Simple & Fast Workflow

Convenient Diagnosis Function

Contact us:

3555 Veterans Memorial Hwy, Unit D
Ronkonkoma, NY 11779
sales@jpihealthcare.com
(516) 513-1330 Ext. 3



www.jpihealthcare.com