



ExamVue DUO



Easy to Learn

Quick to Use

No Hidden Fees

ExamVue Duo Acquisition Software

ExamVue Duo Acquisition Software is JPI Healthcare Solutions's advanced software, designed for the acquisition, processing, storage, and viewing of digitally acquired x-ray images.

ExamVue Duo is indicated for use in general radiology, specialist radiology (podiatry, orthopedic, chiropractic, and other specialties), mobile x-ray applications. Large or small, ExamVue Duo has the tools and ease of use to fit your clinic.

Our Solutions - Your Benefits

Developed and supported in-house by JPI, our total ownership of our software allows us to offer features and support in a way unavailable from most competitors.

JPI is dedicated to making your experience with ExamVue Duo the best it can possibly be. Many features that are expensive "options" in other software are included in every ExamVue Duo install.

Image Stitching, DICOM Modality Worklist, Procedure Code Matching, DICOM Print, and CD Burning are all standard with ExamVue Duo.



New Software Designed for High-end and High Volume Market



Optimized for Simple and Fast Workflow



Includes Tools Valuable for Both Large Companies and Small Offices



Pictorial Instructions and Techniques for All Exposures



Designed for Touchscreen and Desktop




Complete Integration with JPI CubeX Generators plus Leading Manufacturers, such as CPI, Poskom, and IMD

ExamVue Duo

With state of the art intuitive imaging tools, our ExamVue Duo software makes for easy acquisition and short learning curves for acquiring high-quality images on any patient in your work environment.

ExamVue Duo can be used on Windows tablets, PCs, and laptops, and integrates with multiple detectors, industry-leading x-ray generators, PACS, and EMR systems. Reliable and flexible, ExamVue Duo seamlessly integrates into any professional practice.



ExamVue Duo can be expanded with the ExamVue PACS and JPI Cloud PACS for viewing on other computers and adding further tools for medical specialties including Orthopedics, Podiatry, and Chiropractic. The JPI Cloud PACS expands your access to your images to anywhere you can connect to the internet and allows for easy sharing with doctors and patients throughout your practice or anywhere in the world.

Software At A Glance




Image Manipulation and Annotation Tools:


Image Manipulation: Flip (vertical/horizontal), Rotate (clockwise/counter-clockwise), Zoom, Pan, Window/Level, Magnifying Glass, Crop

Image Annotations and Measurements: Line Measurement, Angle Measurement, Free Text, Arrow, L&R Markers, Cobb's Angle Measurement, Norberg Angle Measurement

And many more.

Additional Features:

- ▶ Easy workflow with the fewest necessary interruption during exam acquisition. Complete an entire exam from the first x-ray the last without ever needing to touch the computer.
- ▶ Automatic image processing for all views, plus manual control for fine tuning images.
- ▶ Continual status updates for panel connection, PACS, MWL, and hard drive space.
- ▶ MWL and Procedure Code Matching for integration with RIS/HIS or our stand-alone MWL generator.
- ▶ Export and Import DICOM and JPG files.
- ▶ Integrated technique charts, positioning guides and exposure index for all programmed views.



Changeable Skins



Dual Monitors



Multiple Views

Contact us:

3555 Veterans Memorial Hwy, Unit D
Ronkonkoma, NY 11779
sales@jpihealthcare.com
(516) 513-1330 Ext. 3



www.jpihealthcare.com