

ACQUISITION SOFTWARE

FOR VETERINARY



Easy to Learn

Quick to Use

Specialty Tools Available

ExamVue Duo Acquisition Software

ExamVue Duo Acquisition Software is JPI Healthcare Solutions's advanced X-ray imaging software, designed for the acquisition, processing, storage, and viewing of digitally acquired x-ray images.

ExamVue Duo is indicated for use in veterinary applications. Large or small, ExamVue Duo has the tools for your practice and the ease of use to fit your veterinary clinic.



Developed and supported in-house by JPI, our total ownership of our software allows us to offer features and support in a way unavailable from most competitors.

JPI is dedicated to making your experience with ExamVue Duo the best it can possibly be. Many features that are expensive "options" in other software are included in every ExamVue Duo install.

Image Stitching, DICOM Modality Worklist, Procedure Code Matching, DICOM Print, and CD Burning are all standard with ExamVue Duo.



New Software Designed for the Veterinary Market



Optimized for Simple and Fast Workflow



Includes Tools Valuable for Both Large Companies and Small Offices



Pictorial Instructions and X-ray Techniques for Dogs, Cats, Horses, and Other Animal Exams



Designed for Touchscreen and Desktop



Integration with JPI CubeX Generators plus Leading Manufacturers such as CPI, Poskom, and IMD


ACQUISITION SOFTWARE

FOR VETERINARY

ExamVue Duo

With state of the art intuitive imaging tools, our ExamVue Duo software makes for easy acquisition and short learning curves for acquiring high-quality images on any patient in your work environment.

ExamVue Duo can be used on Windows tablets, PCs, and laptops, and integrates with multiple detectors, industry-leading x-ray generators, PACS, and EMR systems. Reliable and flexible, ExamVue Duo seamlessly integrates into any professional practice.



ExamVue Duo can be expanded with the ExamVue PACS and JPI Cloud PACS for viewing on other computers and adding further tools for veterinary applications. The JPI Cloud PACS expands your access to your images to anywhere you can connect to the internet and allows for easy sharing with doctors and patients throughout your practice or anywhere in the world.

Software At A Glance





Registering

Positioning

An extensive array of tools for measuring and displaying your images, for you to customize to your practice. Tools include, among many others:



Plus optional add-ons for special Companion Animal Tools, Equine tools, or both, with written and graphic instructions on their use.

Additional Features:

- ► Easy workflow with the fewest necessary interruption during exam acquisition. Complete an entire exam from the first x-ray to the last without ever needing to touch the computer.
- Automatic image processing for all views, plus manual control for fine tuning images.
- Continual status updates for panel connection, PACS, MWL, and hard drive space.
- MWL and Procedure Code Matching for integration with RIS/HIS or our stand-alone MWL generator.
- ▶ Export and Import DICOM and JPG files.
- ▶ Integrated technique charts, positioning guides and exposure index for all programmed views.







