



TruePortable

Versatile Imaging. On-the-go.

Lighter, more maneuverable, and better protected than the competition, TruePortable™ was designed specifically with portable provider technologists in mind. The TruePortable™ can be maneuvered easily up stairs, in tight areas, and underneath beds.

Lightweight | Easy-to-use | Seamless Workflow

Our Solutions

With superior image quality and intuitive imaging software, the ExamVue Duo digital imaging system delivers fast image acquisition, easy workflow management, and lower dose for your patients.

The ExamVue Duo system comes equipped with a strong suite of medical tools, and is available in any configuration for your practice needs.

ExamVue Duo Package



ExamVue Duo
14"x17" Wireless



Laptop
Workstation



X-ray Grid



ExamVue Duo
Acquisition Software



DICOM
Modality Worklist

Choose Your PACS



LAN PACS



Cloud PACS

Powerful & Versatile

- ▶ 82lbs/37kg (Generator and Stand)
- ▶ 140µm Csl panel
- ▶ 2.7 kW generator
- ▶ Control techniques directly from laptop

Unmatched Coverage Program

- ▶ Up to 5 years of PANEL WARRANTY
- ▶ Up to 5 years of REMOTE SOFTWARE SUPPORT
- ▶ Up to 5 years of DROP COVERAGE PROGRAM
- ▶ 3 years of LAPTOP WORKSTATION WARRANTY
- ▶ 2 years of X-RAY GENERATOR WARRANTY
- ▶ 1 year of MOBILE STAND WARRANTY






Premium Protection

TruePortable™ combines portability, durability, power, and advanced image processing to provide exceptional images. We built TruePortable™ with protection of your investment in mind. Protective encasements and carry bag protect your detector and laptop during transportation and use.

Great protection will allow your team to go about their day knowing that your most valuable components are well-protected no matter where you take them.



ExamVue Duo Acquisition Software

With state of the art and intuitive imaging tools, our ExamVue Duo acquisition software produces clear, concise, and quality images with even the oldest generators. Reliable and convenient, the ExamVue Duo software seamlessly integrates into any professional practice.



ExamVue PACS comes standard with medical diagnostic tools and can be expanded even further with the JPI Cloud PACS. The JPI Cloud PACS expands your access to your images to anywhere you can connect to the internet and allows for easy sharing with doctors and patients throughout your practice or anywhere in the world.

Ergonomic Design



Folded Vertically

Folded Horizontally

Unfold and be ready to take an X-ray immediately

User-Friendly Design

Simple & Fast Workflow

Convenient Diagnosis Function

Contact us:

3555 Veterans Memorial Hwy, Unit D
Ronkonkoma, NY 11779
sales@jpihealthcare.com
(516) 513-1330 Ext. 3



www.jpihealthcare.com