



TrueMobile APEX

Clarity and Mobility That Transform Bedside Imaging



ULTRA-SHARP IMAGE QUALITY:

99 μm resolution for precise, detailed mobile X-ray images.



LONG-LASTING OPERATION:

Up to 1,200 exposures with a full charge.



8 YEARS OF BATTERY WARRANTY:

Maximum uptime without the need for replacements.

HIGHLIGHTS FEATURES



LIGHTWEIGHT & MANEUVERABLE:

Zip through narrow hallways and tight spaces with ease.



LONG REACH FOR BEDSIDE EXAMS:

Allows proper positioning with limited patient movement.



HOT-SWAP BATTERY:

No downtime during exams.

Fits in most elevators

Unobstructed view, with no blind spots

8-year warranty for the battery
2-year warranty for the system

Lightweight design with powerful performance from a 32 kW generator



TrueMobile Apex delivers rapid performance and fast charging—just 10 minutes of charge provides 1 hour of use, reducing the need for multiple mobile X-ray units to meet demand

ExamVue™ Duo Acquisition Software

ExamVue™ Duo is a powerful X-ray imaging software built to deliver reliable performance and exceptional user support. Designed to streamline every step of the radiographic workflow, ExamVue™ Duo brings speed, precision, and simplicity to modern medical imaging.



Easy to Use and Learn



Patient Positioning Guides



Image Auto-Stitching

Customer Service

We are JPI. **After sales service and support** isn't just a priority—it's the core of our commitment. We are relentless in our pursuit of one singular goal: to **minimize your downtime**. Our outstanding support team of experienced and knowledgeable engineers are always ready to help.

