


**Clarity That Reveals Details
Standard Panels Leave Behind**


Made in France

Trusted Worldwide – **About 50% of all X-ray exams** across the globe are performed with **Pixium panels**, making them one of the **top detectors in the world**.



EXCEPTIONAL IMAGE QUALITY

Trixell Pixium panel with 99 microns.



LONG-LASTING DURABILITY

Glassless, IP67 waterproof, and drop-tested up to 4 ft.



ENHANCED CYBERSECURITY

Designed to help keep data safer, more private, and compliant.


Industry Leading Product Warranty



FAST REPLACEMENT



7 YEAR WARRANTY



7 YEAR DROP INSURANCE

Workflow optimization that keeps your process simple, fast, and efficient

Hot-swap batteries—no need to power down

Wireless and compatible with any X-ray system



661lbs Weight Capacity

Available in 3 sizes:

17" x 17" Wireless
6.6 lbs with battery


14" x 17" Wireless
5.7 lbs with battery

10" x 12" Wireless
3.5 lbs with battery

Long-lasting battery with up to 7 hours of power

ExamVue™ Duo Acquisition Software

ExamVue™ Duo is a powerful X-ray imaging software built to deliver reliable performance and exceptional user support. Designed to streamline every step of the radiographic workflow, ExamVue™ Duo brings speed, precision, and simplicity to modern medical imaging.



Easy to Use and Learn



Patient Positioning Guides



Auto-Stitching



Customer Service

We are JPI. **After sales service and support** isn't just a priority—it's the core of our commitment. We are relentless in our pursuit of one singular goal: to **minimize your downtime**. Our outstanding support team of experienced and knowledgeable engineers are always ready to help.

